

ADVANCED

SPIROLABTM

Desktop, Stand-alone and
PC-based Spirometer,
with Oximetry Option

All-in-one Spirometer with 7" display,
Embedded Printer and Oximetry option,
to carry on the go



MAIN features



REAL-TIME TEST

Spirometry: FVC, VC, IVC, MVV, PRE/POST
Bronchodilator comparison
Oximetry (optional): Spot test (SpO2%, Pulse BPM)



CARRY EVERYWHERE

7" LCD Color Touchscreen Display, Long-lasting rechargeable battery, massive Internal Storage



COMPLIANCE ATS/ERS 2019

And other Standards including ISO 26782 (for Spirometry), ISO 23747 (for PEF), ISO 80601-2-61 (for Oximetry), and more. CE0476, FDA 510 (k)





PC CONNECTION AVAILABLE

Real-time test on PC screen, connect with your EHR/EMR, back-up internal memory and more, via USB and Bluetooth



CALIBRATION SOFTWARE

Available on device, with printable calibration report (no separate software required)



DISTINCTIVE features



PREDICTED SETS & VALUES

Large Selection, including GLI, comparison %Pred, Z-score and LLN



PRINTED EMBEDDED

Thermal printer. Paper size 112mm. Direct external print also available via PC software



PEDIATRIC INCENTIVE

Real-time animation on display to improve patient compliance during the test



COVID-19 PREVENTION

Complete Disposable Set with Antiviral filter. Bluetooth connection to test at safety distance



Always INCLUDED

- ✎ Carrying case
- ✎ Power supply/battery charger
- ✎ USB cable
- ✎ 1 Roll of thermal printer paper
- ✎ Noseclip
- ✎ PC Software license
- With Oximetry Option:
 - ✎ Finger Probe



Compatible SOFTWARE

winspiroPRO



Pediatric Incentive (PATENTED) to improve patient compliance during the test.

Acceptability Messages, Test interpretation and Quality Control Grade according to the latest Spirometry Standards

MAIN FEATURES

Windows-based solution for Spirometry, Oximetry and Telemedicine.


Wide range of predicted sets and values, including GLI Predicted sets, LLN and Z-score.

Embedded EHR/EMR connectivity.

NET VERSION available, share one database between different PC workstations.

MEDICAL REPORT

Specialized and customizable printout



spiro Connect



MAIN FEATURES

Windows-based solution, direct integration with your EHR/EMR.

Real time test include Spirometry and Oximetry Standardized communication in HL7 or Exchange Protocol.

Select patient info directly from your own EHR/EMR

Spirometry test: FVC-Pre, FVC-Post, VC-Pre
Oximetry Test: SpO2 (%), Pulse (BPM)

GO-TO-MARKET TOOLKIT

Software Development Kit available for System Integrators and App Developers.
OEM service available for Spirometry and Oximetry.



Learn more about available SDK and OEM



Compatible TURBINES

flowMIR™ Disposable Turbine



Reusable Turbine



Mouthpiece

Included Disposable

Required, Not Included

Turbine Disinfection

Not required

Required

Turbine Calibration

Not required

Required

Packaging

Individually sealed: 60 or 10 units / box

1 unit in Carton box

Antiviral Filter

Available Disposable

Required Disposable

PLAY VIDEO



SCIENTIFIC PUBLICATIONS



TECHNICAL datasheet

PRODUCT CODES - Spirolab Configurations

911080EO – Spirometer • 911080E1 - Spirometer with reusable turbine

911081EO - Spirometer + Oximeter • 911081E1 - Spirometer + Oximeter with reusable turbine

Technical specification

Width 220 mm
Length 210 mm
Thickness 51 mm
Weight 1450 g (battery pack included)

Sensors



miniflowmeter (code 900595)
for reusable and disposable turbine
dimension (Ø 30 mm, 42 mm)

Reusable soft, adult, MIR sensor for
oximetry tests (code 919024) only
spirolab code 911081

Power supply
power Rechargeable battery and mains
Current capacity Ni-MH, 6 elements
Consumption 4500 mAh
Backup battery voltage average 250 mA
Batteries charger none

Output voltage=12 V, current=1A,
compliant with EN 60601-1
~10 hours

Autonomy
Connectivity USB 2.0, Bluetooth® 2.1
Display 7 inch colour touch screen LCD
Display with 800x480 resolution
Keyboard absent, touchscreen
Mouthpieces Ø 30 mm (1.18 inch)
Type of electrical Internally powered
Protection Class II while charging battery
Safety level for shock hazard Type BF Apparatus
Conditions of use Apparatus for continuous use

Temperature: MIN -40 °C,
MAX +70 °C
Humidity: MIN 10% RH;
MAX 95%RH

Storage conditions

Temperature: MIN -40 °C,
MAX +70 °C
Humidity: MIN 10% RH;
MAX 95%RH

Transport conditions

Temperature: MIN -40 °C,
MAX +70 °C
Humidity: MIN 10% RH;
MAX 95%RH

Operating conditions

Temperature: MIN +10 °C,
MAX +40 °C
Humidity: MIN 10% RH,
MAX 95%RH

Applied norms

Electrical Safety EN 60601-1
Electro Magnetic Compatibility
EN 60601-1-2

Degree of protection
against water penetration

IPX1 appliance protected against
water leaks

Codes and equipments

911080EO spiro
911080E1 spiro with reusable turbine
911080E2 spiro with 120 FlowMir
911081EO spiro+oxy
911081E1 spiro+oxy with reusable turbine
911081E2 spiro+oxy with 120 FlowMir

Spirometry

Flow sensor bi-directional digital turbine
Volume rate 10 L
Flow range ±16L/s
Volume accuracy ±2.5% or 50 mL
Flow accuracy ±5% or 200 mL/s
Dynamic resistance <0.5 cm H₂O/L/s
Temperature sensor semiconductor (0-45°C)
Test available FVC, VC, IVC, MVV, PRE-POST
Measured parameters FVC, FEV1, FEV1/FVC%, FEV1/PEF,
FEV1/VC, FEV1/FEF0.5, DTPEF, FEV
FEV0.5/FVC, FEV0.75, FEV0.75/FVC,
FEV2, FEV2/FVC, FEV3, FEV3/FVC,
FEV1/FEV6, PEF, FEF25, FEF50,
FEF2575, FEF7585, FET, Vext, ELA,
FIVC, FIV1, PIF, FIV1/FIVC, FIF25,
FIF75, R50, MVVcal, PIF, IRV, VC,
IVC, IC, ERV, IRV, FEV1/VC, TV, VE,
te, ti/t-tot, TV/ti, MVV
Memory capacity Up to 10000 tests

Oximetry (on request)

Measurement method Red and infrared absorption
SpO₂ range 0-99%
SpO₂ accuracy ± 2% between 70-99% SpO₂
Average number of heart beats for the %SpO₂ calculation 8 beats
Pulse Rate range 18-300 BPM
Pulse Rate accuracy ± 2BPM or 2% whichever is greater
Average interval for the calculation of cardiac pulse 8 seconds
Signal quality indication 0 - 8 segments on display
Test available spot
Measured parameters SpO₂% min, max, average
BPM min, max, average
Test duration
% Bradycardia Duration (<40 BPM)
% Tachycardia Duration (>120 BPM)
% of Time with SpO₂ ≤ 90% (T90%,
T89%), T5
Memory capacity about 500 hours oximetry

Certificates & Registrations

CE 0476 MED 9826
FDA 510 (k) K 052140
Health Canada 71191 (class II)
CND code Z12150102 (spiro)
Z1203020408 (spiro + oxy)
46906 (spiro), 45607 (spiro + oxy)
GMDN code 1272475/R (spiro)
Ministry of Health 1272476/R (spiro + oxy)
1645455/R (spiro)

ITALY

MIR Head Office
Via del Magliolino, 125
00155 Roma
Tel. +39 06 22 754 777
Fax +39 06 22 754 785
Mir.spirometry.com

USA

MIR USA, Inc.
5462 S. Westridge Drive
New Berlin, WI 53151
Phone +1 (262) 565-6797
Fax +1 (262) 364-2030

FRANCE

MIR Local Office
Jardin des Entreprises,
290, Chemin de Saint Dionisy
30980 LANGLADE (France)
Phone +33 (0)4 66 37 20 68
Fax +33 (0)4 84 25 14 32