



# Spirodoc



**Touchscreen portable spirometer  
with optional 3D oximeter: 6MWT,  
Sleep Test, and 24h Holter for SpO<sub>2</sub>%**

# Available configurations

Spirodoc is available in 3 configurations:

## Spirometer



## Spirometer + Oximeter



## Oximeter



# Supported tests

**Spirometry:** FVC, VC, MVV, PRE/POST bronchodilator comparison

**Oximetry (optional):** Spot test (SpO2%, BPM), 6MWT, Sleep test, and 24h Holter for SpO2%

## Key features

### Touchscreen

Touchscreen for fast data entry

### Multifunctional

In addition to spirometric and oximetric spot tests, Spirodoc also makes it possible to perform 6MWT, Sleep Test, and 24h Holter Test with the 3D oximetry option



### 3D Oximetry: 6MWT, Sleep Test, and 24h Holter for SpO2%


Measurement of desaturation events during exercise, sleep and daily activities

### Optoelectronic reader for removable spirometry

Available in Spirometer and Spirometer + Oximeter configurations to facilitate oximetry testing

### Triaxial accelerometer (with 3D oximeter)


Triaxial accelerometer for recording patient movement and position during measurement



## Real-time tests

Real-time tests displayed on the device and PC screen

## Pediatric incentive



Real-time animation available on **MIR Spiro** software, for improved patient collaboration during the test

## Integrated temperature sensor

Automatic BTPS Conversion

## Long-lasting rechargeable battery

Long-lasting rechargeable lithium battery for extended autonomy in Stand Alone mode



## Predicted values

Wide selection of predicted values including GLI, ERS and others, directly on the device and in PC mode

## EMR/EHR connectivity

Integration via **MIR Spiro** software with EMR/EHR (in HL7, GDT, FHIR, EXCHANGE PROTOCOL)

# Compatible turbines

|                                 |   | Mouthpiece             | Turbine disinfection | Turbine calibration | Packaging                                 | Antiviral filter   |
|---------------------------------|---|------------------------|----------------------|---------------------|---|--------------------|
| FlowMIR® disposable turbine     |  | Disposable included    | Not required         | Not required        | Individually packaged: packs of 60 pieces | Optional           |
| Single patient reusable turbine |  | Required, not included | Required             | Required            | Pack of 1 unit                            | Recommended by ATS |

# How to use

Spirodoc works both in **Stand Alone** mode and connected to the **PC via USB cable**

## MIR Spiro software

- \\ Comprehensive software for spirometry and oximetry
- \\ Designed to be integrated with EMR/EHR
- \\ Complies with the latest ATS/ERS guidelines
- \\ Available for desktop and laptop use
- \\ MacOS and Windows

All MIR professional devices work with **MIR Spiro** software, **the latest generation software** for spirometry and oximetry.



## Platinum Card

To subscribe to a Platinum subscription plan it is necessary to **have the MIR Spiro Platinum Card.**

# Measured parameters

|                                | From MIR Spiro software<br>via connection<br>to the device   | From device<br>in Stand Alone mode   |
|--------------------------------|--|--|
| <b>Spirometry</b>              | FVC, FEV1, PEF, FEF75, FEF25-75, FET, FEV1/FVC, FEV6, FEV1/FEV6, FEF25, FEF50, FIVC, FEV1/VC, ELA, MVV(cal), Time to PEF, FEV0.5, FEV0.5/FVC, FEV0.75, FEV0.75/FVC, FEF75-85, Extr. Vol, VC, EVC, IVC, IC, VC, ERV FEV3, FIV1, FIV1/FIVC, PIF, FEV3/FVC, PIF, FEV2, FEV2/FVC, FIF25, FIF50, FIF75, R50, FEV1/PEF (EI), FEV1/FEV0.5 (RFEV), TV, VE, RR, tI  | *FVC, *FEV1, *PEF, FVC, FEV1, FEV1/FVC, FEV1/VC, PEF, FEF25-75, FEF25, FEF50, FEF75, FEV3, FEV3/FVC, FEV6, FEV6%, FET, BEV, FIVC, FIV1, FIV1/FIVC, PIF, MVVcal, VC, EVC, IVC, IC, ERV, TV, VE, RR, tI, tE, TV/tI, tI/tTOT, MVV, ELA<br>*Best values  |
| <b>Oximetry<br/>(optional)</b> | SpO2% [Baseline, Min, Max, Mean], BPM [Baseline, Min, Max, Mean], T Total, T Analysis, T <90%, T <89%, T <88%, T <87%, EvSpO2% <89, Δ Index, t <40BPM, t >120BPM, Ev <40BPM, Ev >120BPM, SpO2% End, BPM End, SpO2% Start, BPM Start, T Walk, T Recovery, Distance, T2%Δ SPO2, T4%Δ SPO2, Theoretical, Theoretical min, Theoretical %, Theoretical min, AUC/Distance*, Dyspnea Start, Dyspnea Fin, Dyspnea CHG, Diastolic Start, Systolic Start, Diastolic End, Systolic End, Steps, VMU**, O2-GAP***, O2, ODI Average Desat., Tot Desatur., Max Duration, Peak Desatur., BPM Index, Average Desat., Average Fall, Max Fall, BPM Change, NOD4%, NOD89%, NOD90%, t.NOD4%, t.NOD89%, t.NOD90% | SpO2% [Baseline, Min, Max, Mean], BPM [Baseline, Min, Max, Mean], T Total, T Analysis, T <90%, T <89%, T <88%, T <87%, EvSpO2% <89, Δ Index, t <40BPM, t >120BPM, Ev <40BPM, Ev >120BPM, SpO2% End, BPM End, SpO2% Start, BPM Start, T Walk, T Recovery, Distance, T2%Δ SPO2, T4%Δ SPO2, Theoretical, Theoretical min, Theoretical %, Theoretical min, AUC/Distance*, Dyspnea Start, Dyspnea Fin, Dyspnea CHG, Diastolic Start, Systolic Start, Diastolic End, Systolic End, Steps, VMU**, O2-GAP***, O2, ODI Average Desat., Tot Desatur., Max Duration, Peak Desatur., BPM Index, Average Desat., Average Fall, Max Fall, BPM Change, NOD4%, NOD89%, NOD90%, t.NOD4%, t.NOD89%, t.NOD90% |

# Datasheet

code 911080xx (spiro) code 911081xx (spiro+oxy)

|                                      |  |
|--------------------------------------|--|
| <b>Main body</b>                     |  |
| <b>Size</b>                          | 48 x 101 x 16 mm   |
| <b>Weight</b>                        | 99 g (battery included)  |
| <b>Turbine housing</b>               |  |
| <b>Size</b>                          | 47 x 46 x 24 mm  |
| <b>Weight</b>                        | 17 g (battery included)  |
| <b>Turbines</b>                      | <ul style="list-style-type: none"> <li>· Reusable turbine (code 910002)</li> <li>· Disposable turbine (code 910004)</li> </ul>   |
| <b>Accelerometer</b>                 | triaxial accelerometer   |
| <b>Power supply</b>                  | 3.7V lithium-ion battery, 1100 mAh rechargeable  |
| <b>Current</b>                       | 1100 mAh   |
| <b>Consumption</b>                   | -20-30 mA (during testing)   |
| <b>Charge Batteries</b>              | Voltage =5 V DC,<br>Current = minimum 500 mA,<br>Connector: micro USB type B<br>Complies with EN 60601-1   |
| <b>Autonomy</b>                      | 50 hours   |
| <b>Connectivity</b>                  | USB 2.0, Bluetooth® 2.1  |
| <b>Display</b>                       | monochrome LCD,<br>160 x 80 pixels<br>Size 2.8 inches  |
| <b>Keyboard</b>                      | Touchscreen  |
| <b>Mouthpiece</b>                    | Ø 30 mm (1.18 inches)  |
| <b>Type of electrical protection</b> | Powered internally   |
| <b>Safety level</b>                  | Type BF device   |
| <b>Against shock</b>                 |  |
| <b>Terms of use</b>                  | Device for continuous use  |
| <b>Storage conditions</b>            | Temp: MIN -20°C, MAX+60°C<br>Humidity: MIN 10% RH; MAX 95%RH   |
| <b>Operating conditions</b>          | Temp: MIN +10°C, MAX +40°C<br>Humidity: MIN 10% RH, MAX 95%RH  |
| <b>Applicable regulations</b>        | Electrical Safety IEC 60601-1<br>Electro Magnetic Compatibility EN 60601-1-2<br>ISO 80601-2-61:2017<br>ISO 26782: 2009<br>ISO 23747: 2015<br>ATS/ERS:2005, 2019(update)<br>IEC 60601-1-6:2010<br>IEC 60601-1-8:2006+ AMD1:2012<br>IEC 60601-1-9:2007+AMD1:2013<br>IEC 62304:2006 + A1:2015<br>ISO 10993-1:2018<br>Directive 2014/53/EU RED |

|  |  |
|--|--|
| <b>Spirometry</b>                                    |  |
| <b>Sensor</b>  | two-way digital turbine  |
| <b>Flow range</b>                                    | ±16L/s   |
| <b>Volume accuracy</b>                               | ±2.5% or 50mL  |
| <b>Flow accuracy</b>                                 | ±5% or 200 mL/s  |
| <b>Dynamic resistance</b>                            | <0.5 cm H2O/L/s  |
| <b>Temperature sensor</b>                            | semiconductor (0-45°C)   |
| <b>Measured parameters</b>                           | FVC, FEV1, FEV1/FVC%, FEV3, FEV3/FVC%, FEV6, FEV1/FEV6%, PEF, FEF25%, FEF50%, FEF75%, FEF25-75%, FET, Vext, ELA, FIVC, FIV1, FIV1/FIVC%, PIF, VC, IVC, EVC, IC, ERV, FEV1/VC%, VT, VE, RR, ti, te, ti/t-tot, VT/ti, MVV                |
| <b>Memory capacity</b>                               | more than 10,000 tests   |
| <b>Oximetry (on request)</b>                         |  |
| <b>Measurement method</b>                            | Infrared absorption  |
| <b>SpO2% Range</b>                                   | 0-99%  |
| <b>Accuracy of SpO2%</b>                             | ± 2% between 70-99% SpO2   |
| <b>Average number of beats for SpO2% calculation</b> | 8 beats  |
| <b>Cardiac pulse range</b>                           | 30-254 BPM   |
| <b>Cardiac pulse accuracy</b>                        | ± 2BPM or 2% the greater of the two  |
| <b>Mean interval for calculation of heartbeat</b>    | 8 seconds  |
| <b>Signal quality indication</b>                     | 0 - 8 segments on screen   |
| <b>Measured parameters</b>                           |  |
| <b>For each test</b>                                 | SpO2%MIN, SpO2%MEAN, SpO2%MAX, BPM MIN, BPMMEAN, BPM MAX, Ttotal, Tanalysis, T<90%, T<89%, T<88%, T<87%, EvSpO2%<89, ΔIndex, T<40BPM, T>120BPM, Ev<40BPM, Ev>120BPM  |
| <b>Sleep Test</b>                                    | SpO2%BASE, BPMBASE, ODI, Mean Dur. Desat., TotDesaturat., Longest Desat., Desatur. Peak, BPM Index, Mean Desaturat., Mean Drop, Max Drop, BPM Variation, NOD4%, NOD89%, NOD90%, t.NOD4%, t.NOD89%, t.NOD90%<br>Record of body position |

|                                       |   |
|---------------------------------------|---|
| <b>6MWT Test</b>                      | SpO2% start, SpO2% end, BPM start,BPM end, SpO2% base, Tbaseline, Twalking, Trecovery, Distance, T2%ΔSPO2, T4%ΔSPO2, Predicted, %Predicted, Predicted min, %Predicted min AUC/Distance, Dyspnea, Dyspnea base, Dyspnea end, Dyspnea CHG, Fatigue base, Fatigue end, Fatigue CHG, Diastolic base, diastolic end, Systolic base, Systolic end, Steps, VMU, O2-GAP, O2 |
| <b>Memory capacity</b>                | about 300 hours of oximetry   |
| <b>Certificates and registrations</b> |   |
| <b>CE 0476</b>                        | MDR 2017/745  |
| <b>FDA 510 (k)</b>                    | K 103530  |
| <b>Health Canada</b>                  | 71191 (class II), 75535 (class III)   |
| <b>EMDN liv.4</b>                     | Z121501   |
| <b>CND Code</b>                       | Z12150102 (spiral)<br>Z1203020408 (spiro + oxy)   |
| <b>GMDN Code</b>                      | 46906 (spiral),<br>45607 (spiro + oxy)  |
| <b>Ministry of Health</b>             | 2493989/R (910600)<br>2494292/R (910606)<br>2494301/R (910610)<br>2494198/R (910600I1)<br>2494295/R (910606I)<br>2494319/R (910610I1)<br>2494380/R (910600I0)<br>2494386/R (910610I0)   |

## Compliance with guidelines and standards

**Spirometry:** ATS/ERS 2005 + update to 2019;

ISO 23747: 2015; ISO 26782: 2009

**Oximetry:** ISO 80601-2-61:2017



## ITALY

MIR Medical  
International Research  
S.p.A.

Viale Luigi Schiavonetti,  
270 00173, Rome

Tel. +39 06 22 754 777

Fax +39 06 22 754 785

[mir@spirometry.com](mailto:mir@spirometry.com)

[spirometry.com](http://spirometry.com)

## USA

MIR USA, Inc.  
5462 S. Westridge Drive  
New Berlin, WI 53151

Tel. +1 (262) 565-6797

Fax +1 (262) 364-2030

[mirusa@spirometry.com](mailto:mirusa@spirometry.com)

## FRANCE

MIR Local Branch  
Jardin des Entreprises, 290,  
Chemin de Saint Dionisy  
30980 LANGLADE

Tel. +33 (0)4 66 37 20 68

Fax +33 (0)4 84 25 14 32

[mirfrance@spirometry.com](mailto:mirfrance@spirometry.com)

## BRAZIL

MIR Local Branch  
Rua Pinheiro Machado, 2659,  
Sl.303, Caxias do Sul RS

Tel +55 5430253070

[mirbrazil@spirometry.com](mailto:mirbrazil@spirometry.com)

