

Winspiro PRO™ & Winspiro PRO NET™

Complete and powerful software package for Specialists.
Includes Telemedicine and Home Care Management

Winspiro PRO™

Customizable PC software for Spirometry,
Oximetry and Telemedicine projects

*FVC, VC, IVC, MVV, Bronchodilator Reversibility, Bronchial
Challenge with FEV1-response curve with protocols for both
Methacholine and Mannitol*

Paediatric Incentive for Spirometry

Home Care and Telemedicine management

Quanjer 2012-GLI Predicted values with LLN and Z-score

Patient Trend Charts for easy follow-up

Search engine for instant access to data

Specialized and customizable Printout

Wide choice of communication protocols ... and more

3D Oximeter® management, O2Gap Index,

Sleep Desaturation Analysis, 6 Minute Walk Test, Daytime

Physical Activity Report coupled with Desaturation Events


Winspiro PRO NET™

Network version of Winspiro PRO™


PC software for Spirometry,
Oximetry and Telemedicine

*Complements the features of the client version,
by providing the capability to share a single
Database between different network users*

Winspiro PRO™ is compatible with Windows: XP, VISTA, 7, 8 and 10.
Winspiro PRO™ includes free updates and is provided with Spirolab™,
Spirodac™, Spirobank II™, SpiroTel™, Minispir™, Spirolab III™



Screenshot summary of all tests carried out



6 Minute Walk Test:
Events and Physical Activity recording

Winspiro PRO NET™ includes free updates, and is provided on request