

Spirometry Results




VISIT DATE 2/2/2021

ID	1234567890	Origin	Caucasian	Weight	210 lb
Last Name	Test Patient	Age	66	BMI	30.1
First Name	Sample	Gender	M		
Date of birth	3/13/1957	Height	5' 10		

ACCEPTABILITY CRITERIA

Quality Grade PRE A Variability FEV1=0.01(0.43%), FVC=0.22(7.83%)

Acceptable trials 2



Spirometry

Parameters		Best	LLN	Pred.	%Pred.	PRE#1	PRE#2	PRE#3
FVC	L	3.03	3.27	4.37	69.3	3.03	2.81	2.68
FEV1	L	2.32	2.44	3.33	69.6	2.31	2.32	2.24
FEV1/FVC	%	76.57	63.51	76.45	100.2	76.2	82.6	83.6
PEF	L/m	588	300	504	117.2	558	540	588
FEF2575	L/s	1.83	1.18	2.61	70.2	1.83	2.39	2.41
FEF25	L/s	5.16	0.0	0.0	-	5.16	5.92	6.7
FEF50	L/s	2.1	0.0	0.0	-	2.1	2.49	2.5
FEF75	L/s	0.42	0.28	0.74	57.1	0.42	0.71	0.87
PEFTime	ms	24	-	-	-	24	30	19
EVol	mL	44.0	-	-	-	44.0	52.0	46.0
FEV6	L	2.99	3.27	4.37	68.4	2.99	2.81	2.68

BTPS 1.16 0.0 °C (32.0 °F) - Pressure: 0.0 atm

Oximetry

SpO2 mean 95.0% (min 93.0% - max 96.0%) Pulse mean 95.0 BPM (min 93.0 - max 97.0)

NOTE

Sample test performed on the MIR SBS with Oxi.